Discover the easy way to hear better with digital technology









How it works

Maxi's built-in microphone picks up and amplifies sounds using digital technology. You hear the sounds through headphones, earphones or a neck loop.

Maxi also has a built-in telecoil that works with inductive loop systems that are common in public areas.





Headphones - for easy use

Stereo headphones with superior sound quality. The headphones are easy to use and have a robust and attractive design. *Model: BE9122*



Earphones - for a discrete look

Stereo earphones with excellent sound. Light comfortable earphones that are easy to use.

Model: BE9124



Neck loop - for your hearing aid

The neck loop sends sound to hearing aids equipped with a telecoil. Put Maxi close to a speaker or connect it to the TV with the cable kit. Maxi works as a remote amplified microphone and makes the sound clearer.

Model: 8E9128



Cable Kit - for HiFi TV sound

Use the cable kit to connect your Maxi to the TV and enjoy high quality sound in full stereo.

If you push the M-button while listening, you can also hear people talking without removing your listening device.

Model: BE9126



External Microphone

For meetings and easy TV-listening
Specifically developed for Maxi as an easy-to-setup
hearing solution in mono. Place the microphone
close to the TV's loudspeaker. The external microphone can also be used at dinners and conferences to improve hearing.

Model: 8E9127



Bellman Audio Maxi

Whether for your own use or as a gift, Maxi comes with everything you need in a hand-some box. Choose Maxi as your solution and you will discover how easily you can improve your hearing.

Maxi comes with

- Maxi Digital Communication Aid
- Listening device of your choice (headphones, earphones or neck loop)
- Audio cable kit for TV
- Two-year warranty

Contact us

Contact your authorized Bellman Audio Maxi dealer (see reverse side) or Bellman & Symfon® for information about a retailer close to you.

This is your solution



Swedish hearing solutions

Bellman Audio series comprises solutions For 20 years, Bellman & Symfon of Sweden unparalleled ease of use of our products. ologists, users and leading design experts. developed in close collaboration with audithat have been thoroughly researched and life for those with hearing deficiencies. The has worked hard to improve the quality of This guarantees great sound quality and

for everyday life

Bellman & Symfon"

Superior sound

Discover Bellman[®] Audio Maxi

a result, Maxi has crystal clear sound and can be used for hours without you growing tired of it. nology as found in modern CD/DVD players. As cated handling. The sound technology in Maxi, Do you remember the old cassette tapes? Along with the music came unwanted noise and complihowever, is fully digital and uses the same tech-

Easy to use

paired vision or dexterity problems. trols are simple to use, even for those with imberized material makes it easy to grip. The connicely on the table or in your hand, and the rub-Maxi is exceptionally user friendly. Its shape fits

and comes with a two-year warranty. Maxi is designed by Bellman & Symfon of Sweden

Head Office ellman & Symfon Europe AB

China Office

ax: +46 31 68 28 90 hone:+46:31:68:28:20

21 32 Vastra Frolunda

39 Huanshi Dong Rd

angzhou, China, 510098

108, Main Tower, GD Int'l Hote

llman & Symfon Asia Ltd.

Långebergsgatan 30

www.bellman.com -mail info@bellman.com

ax: +86 20 8333 5528 -mail info@bellman.com

© 2006, Bellman & Symfon AB, Bellman", Bellman" Visit and Bellman & Symfon" by Bellman & Symfon" AB, all rights reserved. All Product designs are paterized.

Bellmanaudio

Are

these situations familiar?

The easy way to hear better



Having difficulty hearing?



When talking one-on-one?

when chatting with a friend or relative? Do you find it hard to follow a conversation



at family dinners or get-togethers? Do you have trouble hearing what is being said





mai volume? Do you miss words when you watch TV at nor-



ment, for instance, in your car? Do you find it hard to hear in a noisy environ-



